



Explore



Analyse



Discovery



Models



Review



Library



Jobs · 3



Settings

PROJECT · recorded with runs

Data

Datasets

Channels & events

Vocabulary

Analysis

Nulls

Analysis defaults

Compute & HPC

Blocks

Workspaces

Review queues

Models & registration

Library groupings

Files

Storage & backups

Export

Record

Audit log

About

PERSONAL · this browser

Interface

Display

Keyboard & behaviour

differs from default

Import a recording

dry run · nothing written yet



1 File read, not copied · stays where it is

path

D:/recordings/M5_sep_concat_fs1.mat

format

MATLAB v7.3 · 16 variables

2 Sampling rate read from the file header when present; entered rates are marked inferred

fs

1 Hz

read

samples per channel

1,080,000

duration

300 h

3 Channel map 16 variables → 16 channels · names follow CHn_XY

variable

channel

electrode · position

shared ground

include

ch01

CH1_A1

A1 tip · 0 mm

–

☒

ch02

CH2_A1

A1 base · 12 mm

–

☒

ch03

CH3_A2

A2 tip · 0 mm

CH4_A2

☒

... 13 more

4 Start time and link a linked recording is the same signal at another rate

start

2025-09-01 10:00

time zone

Europe/London

link to existing

none

5 Checks must pass before Import

- ✓ samples = duration × fs for every channel (1,080,000 = 300 h × 1 Hz)
- ✓ channel names unique and not already used by this recording
- ✓ not the held-out recording · M4_aug is locked in this page
- ⚠ species not set · allowed, marks the recording provisional

will create 1 recording · 16 channels · 16 .npy files · 8.2 GB on disk · no runs affected

Cancel

Import

Reset page to defaults

Import a recording...

linked

status

fs2

in use

fs1

in use

–

in use

–

provisional

–

held out · locked

electrode config

sub-dermal pairs · 12 mm

noise floor ⓘ

0.10

mV

notes

timed events → Channels & events

on · M4_aug held out

it · turning off needs the name typed · logged